


FUNCTION

MS Medical Ceiling Pendant Function Distribution of medical gases and electricity in Operation theatres for surgeons and anesthetists.

The Pendant is designed to carry up to 30k.G of equipment-electro-surgical Unit, monitor, etc.

CONSTRUCTIONS

- 2distribution ducts in anodized extruded aluminum length 0.60m
- Horizontal plane rotation adjustable up to 340°
- Optional accessories on Request including:
 - Medical gas outlet (O₂ - N₂o - Air - Vacuum - Surgical Air)
 - Electrical Socket 220V
 - Option of scavenging outlet (AGSS type)



90901915



Can Be Supplied
Upon any Special
Requirements

90901914