

The purpose of MS Medical Humidifier is to humidify the oxygen administered to the patient.

It is then connected to an oxygen mask:

- Unbreakable Polycarbonate jar
- The Humidifier is available in 250cc
- Autoclavable at 121°C & 134°C or cold sterilizable
- Humidifier has a silicone tubing for the diffuser 14cm long



90901904



90901905

New MS Humidifier:

- Unbreakable Polycarbonate jar
- Available in 250cc
- Autoclavable at 121°C & 134°C or cold sterilizable

HUMIDIFIER: TECHNICAL VIEW

1	Spare 250cc bottle
2	Diffuser
3	Silicone tubing for diffuser (14cm long)
4	Nut to tubing
5	O-ring for humidifier
6	Inlet O-ring for Humidifier
7	White Cover
8	Inlet Connector
9	Nut
10	Inlet gasket
11	The mask Connector

