


# Medical Bed Head Unit

MS Medical Bed head unit is built around a Specially profiled aluminum Section which provides Linear rigidity for mounting on imperfect walls and intensive use.



90901913



- Unless otherwise required, the unit are Supplied in standard Lengths of 1.50, 3.00 and 4.5m for Patient rooms with one , Two or Three beds.
- Our bed head unit has individual compartments so that the gases are separate from the electric services.
- Can be quickly installed without being dismantled.
- Our bed head unit available in several colors.
- There are Two lights in MS medical Bed Head unit:
  - 1- Direct Light Lamp (60 cm 18 w)
  - 2- indirect light lamp (120cm 36w)
- The Bed Head unit accessories:
  - Direct Light switch.
  - Indirect Light Switch.
  - Electricity Sockets.
  - Telephone Socket.
  - Nurse Call.
  - Medical gas outlets.

Can Be Supplied  
Upon any Special  
Requirements