

Ms Medical Alarm Device for gas control for the signalization of failure in gas supply and the recharge over on emergency reserves.



90901922

In case of failure in gas the alarm unit will make light and voice alarm. There are a lot of light kind:

Green	Light on-work
Red	Light Stop
Cut red light	Alarm Choice

- Switch for stop the alarm.
- The electrical power of the alarm unit 140 V-50 Hertz.
- There are six models from alarm device.

90901922/1	ALARM UNIT (1 LINE)
90901922/2	Alarm unit (2 line)
90901922/3	Alarm unit (3 line)
90901922/4	Alarm unit (4 line)
90901922/8	Alarm unit (8 line)
90901922/16	Alarm unit (16 line)

