

MS Medical Vacuum jars are used to collect liquids or mucus from the medical suction.





Our vacuum jars are all unbreakable and Autoclavable (121°C & 134°C) with capacity of 2 L.

All Our vacuum jars are fitted with an extremely simple Anti over flow safety device. Our 2L Polycarbonate Jars have a moulded hanger making the use of a basket unnecessary.

- Our vacuum jars are Autoclavable 121°C & 134°C.

- Our vacuum jars have the advantages of easy use and MS medical offers you wide range of accessories:

Metal Baskets, floor stands, universal rail adaptors and Aluminum wall brackets.



**NEW MS MEDICAL VACUUM JAR**

