

MS Medical Bed head unit is built around a Specially profiled aluminum Section which provides Linear rigidity for mounting on imperfect walls and intensive use.



- Unless otherwise required, the unit are Supplied in standard Lengths of 1.50, 3.00 and 4.5m for Patient rooms with one , Two or Three beds .
- Our bed head unit has individual compartments so that the gases are separate from the electric services.
- Can be quickly installed without being dismantled.
- Our bed head unit available in several colors.
- There are Two lights in MS medical Bed Head unit
 - 1-Direct Light Lamp (60 cm 18 w).
 - 2- in Direct light lamp (120cm 36w)
- The Bed Head unit accessories:
 - Direct Light switch .
 - In Direct Light Switch .
 - Electricity Sockets.
 - Telephone Socket.
 - Nurse Call.
 - Medical gas outlets.