


Ceiling Pendant (110)



110-A



110-B



110-C

Function

MS Medical Ceiling Pendant Function Distribution of medical gases and electricity in Operation theatres for surgeons and anesthetists. The Pendant is designed to carry up to 30k.G of equipment-electro-surgical Unit, monitor, etc.

Constructions

- 2distribution ducts in anodized extruded aluminum length 0.60m.
- Horizontal plane rotation adjustable up to 340°
- Optional accessories on Request including:
 - medical gas outlet (O₂-N₂o -Air - Vacuum - Surgical Air)
 - Electrical Socket 220V
 - Option of scavenging outlet (AGSS type)